


# ELVeS<sup>®</sup> smart

## Pullback device



Continuous and homogenous withdrawal of the ELVeS<sup>®</sup> Radial fiber family during the endovenous laser ablation of varicose veins.

ELVeS<sup>®</sup> smart, your add-on for a perfect match using ELVeS<sup>®</sup> Radial fibers and our reliable biolitec<sup>®</sup> diode lasers with ELVeS<sup>®</sup> Signal Mode.



## ELVeS<sup>®</sup> smart

REF	L-ELVesSmart-Set
Pulling speed	0.5; 1.0; 1.5; 2.0 mm/s
Fibers	ELVeS Radial <sup>®</sup> , ELVeS Radial 2ring <sup>®</sup> ELVeS Radial <sup>®</sup> slim, ELVeS Radial <sup>®</sup> 2ring slim
Power Supply	100 - 240 VAC, 50 - 60 Hz
Dimensions (H x W x D)	12 cm x 16,5 cm x 22 cm
Weight	2,4 kg